

# biolutoxR: An R-Shiny package for easy performing data analysis of a toxicity test based on bacterial bioluminescence inhibition

Le Picard Coralie<sup>1,2,3</sup>\*, Cachot Jérôme<sup>2</sup>, Clérandeau Christelle<sup>2</sup>, Bringer Arno<sup>3</sup>, Bellier Benjamin<sup>3,4</sup>

<sup>1</sup>LIENSs, UMR 7266, La Rochelle Université, CNRS, 2 rue Olympe de Gouges, F-17042 La Rochelle Cedex 01, Nouvelle-Aquitaine, France

<sup>2</sup>EPOC UMR 5805, Université de Bordeaux, CNRS, Bordeaux INP, 33600 Pessac Cedex, Nouvelle-Aquitaine, France

<sup>3</sup>Qualyse, syndicat mixte, 5 All. de l'Océan, 17000 La Rochelle Cedex, Nouvelle-Aquitaine, France

<sup>4</sup>INRAE Nouvelle-Aquitaine Bordeaux Centre, UR EABX, 50 Avenue de Verdun, 33612, Cestas Cedex, Nouvelle-Aquitaine, France

## Abstract.

In the 21st century, in the context of an environmental crisis, chemical pollution has become a major global concern. In addition to chemical analysis, many bioassays have been developed and have proved to be interesting tools for considerably improve our understanding of the effects of human activities on species and ecosystems. Compared with animal or plant bioassays, bacterial bioassays based on bioluminescence inhibition have the advantage of being relatively economical and reproducible. This toxicity test uses the bioluminescent metabolic response of bacteria exposed for a short time to a solution of interest in the aim of assessing and quantifying its toxicity. Although routine tools are available to pre-process the data obtained, to our knowledge no flexible research tool is available for the data analysis following these tests, which could limit the access to the results for novices or experienced people unfamiliar with data processing software. To overcome this lack of tools for this bioassay, an R-Shiny package is proposed to generalise data analysis following a toxicity test based on bacterial bioluminescence inhibition. The traditional paper-based working environment is reproduced digitally in this package, which ultimately facilitates data entry and cleaning, makes the creation of relevant dynamic graphs, and simplifies access to toxicity data (e.g. dose-response curve and median effective concentration, i.e. EC<sub>50</sub>). The aim of this tool is to provide the target community with a high-performance tool that can be used to obtain toxicity test results based on the inhibition of bacterial bioluminescence.

## Code metadata.

Current Code Version	version 0.1.0
Permanent link to code / repository used for this code version	<a href="https://github.com/bbellier/biolutoxR_package">https://github.com/bbellier/biolutoxR_package</a>
Permanent link to reproducible capsule	<i>NA</i>
Legal code license	GPL ≥ 3
Code versioning system used	Git
Software code languages, tools and services used	R
Compilation requirements, operating environments and dependencies	R (≥ 4.1.2); dependencies: bslib, DT, openxlsx, shiny, shinythemes, shinyMatrix, tidyverse
If available, link to developer documentation/manual	<a href="https://bbellier.github.io/biolutoxR_website/">https://bbellier.github.io/biolutoxR_website/</a>
Support email for questions	coralie.lepicard@yahoo.com

**Keywords.** Environmental pollution, Chemicals, Ecotoxicology, Bioassay

## 1. Motivation and significance

Six of the nine planetary limits - limits that cannot be exceeded to ensure the processes necessary to preserve the stability and resilience of the earth system - have been exceeded, including the limit for chemicals [1]. Indeed, for more than two centuries, human domestic and industrial activities have contributed to an increase in the discharge of chemical substances into aquatic environments, contributing to a deterioration in the functioning of ecosystems [2,3]. In response to the ecological emergency caused by environmental pollution, standardised methods have become essential for identifying toxic effects and, ultimately, helping to ensure the resilience of the environment [4]. Among the existing standardised methods, chemical analysis of aquatic environments appeared to be one of the most relevant, but ultimately proved to be expensive and limited, which finally complicated routine monitoring [5–7]. Thus, the use of bioassays - experimental methods used to study the toxic effect of a substance on model organisms - has proved to be complementary and less restrictive [8].

Nowadays, numerous bioassays are used, using different types of organisms: animals (e.g. mammals, fish, invertebrates), photosynthetic organisms (e.g. plants, algae) and microbiological organisms (e.g. bacteria, fungi). However, bioassays performed on animal and photosynthetic species have been confronted to some major disadvantages (costly experiments, large sample volume, long manipulation time, standardisation of the species, specific equipment, qualified operator and lack of reproducibility), which is reduced for bioassays on bacteria, offering the possibility of assessing the toxicity of a solution of interest with speed, simplicity and reproducibility [9].

Bacterial bioassays using bioluminescence inhibition assess the bioluminescent response of bacteria following exposure to a solution of interest at a given exposure time, generally a short time, i.e. 5, 15 or 30 minutes [10]. This type of test uses particularly some bacteria, among others : *Aliivibrio fischeri* (e.g. Backhaus and Grimme [11]), *Photorhabdus luminescens* (e.g. Masner et al., 2017 [8]), *Photobacterium leiognathi* (e.g. Muneeswaran et al., 2021 [12]) and *Vibrio harveyi* (e.g. Thomulka et al., 1997 [13]). Under optimum conditions, these bacteria produce bioluminescence, i.e. an enzymatic reaction catalysed by luciferase which releases light (see Supplementary materials n°1).

This bioluminescence is directly linked to the respiratory metabolic process of bacteria [14]. In other words, in a favourable environment, these bacteria emit light. However, when a bacterium is exposed to a toxic substance, this metabolic process is disrupted, resulting in the inhibition of bioluminescence. This reaction is therefore used in this type of test to study the direct link between light intensity and the level of toxicity of a sample of interest compared with a control sample. This inhibition is often calculated as follows (e.g. Si et al., 2024 [15]; Westlund et al., 2018 [16]) :

$$I_{\text{biolu}} (\%) = (B_{\text{Ct}} - B_{\text{St}}) / B_{\text{Ct}} * 100$$

Where  $I_{\text{biolu}}$  represents the inhibition of bioluminescence of the bacteria studied (in %);  $t$  represents the exposure time (in min);  $B_{\text{Ct,rb,t}}$  represents the average of the bioluminescence values of the bacteria exposed, during a time  $t$ , to the control solution; and  $B_{\text{St}}$  represents the average bioluminescence values of the bacteria exposed, during a time  $t$ , to the solution of interest.

76 The toxicity of the solution of interest is then visualised on a dose-response curve, providing a visual  
77 representation of the percentage of inhibition (in %) as a function of the log dilution (in %). This curve can  
78 be used to obtain a reference toxicity value: the  $EC_x$ , which is the effective concentration of a toxicant  
79 causing an X% reduction in light intensity of bacteria [17]. The equation for the 5-parameter dose-response  
80 curve is as follows :

81

82 •  $y = c + (d-c) / (1 + \exp(b * \log(x) - \log(e)))$

83

84 Where y represents the inhibition of bioluminescence of the bacteria studied (in %); x represents dilution  
85 (in %); b represents the curvature parameter; c represents the first plate of the curve; d represents the last  
86 plate of the curve; e represents the  $EC_{50}$ .


87

88 More specifically, and in a non-exhaustive way, the toxicity test based on the bacterial bioluminescence  
89 inhibition can be performed experimentally as follows (Fig. 1): (i) creation of a matrix document to record  
90 the position of a sample on the bowl holder, for subsequent data retrieval; (ii) sample preparation: negative  
91 control, positive control, blank and substance of interest (see Table 1 for an explanation of the samples);  
92 (iii) rehydration and acclimatisation of the lyophilized bacteria; (iv) exposure of the bacteria to each  
93 sample; (v) bioluminescence measurements with a luminometer for each sample and for each exposure time  
94 (generally: 5 min, 15 min and 30 min); and (vi) fill in the matrix document, created previously, with the  
95 bioluminescence values. For statistical reasons, the experiment should include several ‘biological  
96 replicates’ for each sample of the substance of interest (at least 3) and several ‘sets’ (i.e. the whole  
97 experiment should be repeated several times). Many studies use this type of approach, resulting in the  
98 development of commercial solutions. (e.g. Microtox®, CheckLight, BioTox, see Kokkali & Van Delft,  
99 2014 [18]).

100

101 The toxicity test based on the bacterial bioluminescence inhibition therefore produces a large amount of  
102 data in a very short time, with a high level of reproducibility. However, at the end of the test, a significant  
103 number of bioluminescence values need to be recorded and processed to obtain the final test results. Some  
104 machines can be used to produce pre-processed results, sufficient for repeated routine analysis, but offering  
105 little freedom to the user, particularly, for research work. To our knowledge, there is currently no package  
106 or application in the R environment that facilitates the analysis of data following this test. For these reasons,  
107 a Shiny-R package was created, simplifying data entry and automating data processing in an application  
108 with an interface specially designed for toxicity tests based on the bacterial bioluminescence inhibition.  
109 Using the conventional annotation format, this Shiny-R package facilitates data entry, automates data  
110 cleaning and makes results visualisation dynamic, particularly for obtaining the dose-response curve and  
111 reference toxicity values (i.e.  $EC_x$ ), which are often complex and subject to errors. In the long term, this  
112 tool could be generalized/standardised for the data analysis of toxicity tests data based on bacterial  
113 bioluminescence inhibition, with the aim of increasing its reproducibility by its widespread use.

114



115  
 116 **Figure 1.** Main steps toxicity tests based on bacteria bioluminescence inhibition. The negative control is  
 117 represented by a '-', the positive control is represented by a '+', the blank is represented by a 'b' and the  
 118 solution of interest is represented by an 's'. The steps represented by letters in brackets (i, ii, iii, iv, v, vi)  
 119 are explained in the 'Motivation and significance' section.

120  
 121  
 122 **Table 1.** Description of samples used in toxicity tests based on bacteria bioluminescence inhibition.

Dependencies (code used in Fig.1)	Description
Negative control (-)	Non-toxic sample used as a reference to show normal bacterial behaviour.
Positive control (+)	Toxic sample used to confirm correct detection and determination of toxicity by the test.
Blank (b)	Non-toxic solution that has been subjected to the same operations as the samples of the solution of interest. It provides assurance that the steps performed during manipulation do not introduce toxic biases. It can be diluted to assess its dose-response relationship.
Solution of interest (s)	Solution of interest, i.e. an ecologically/ecotoxicologically relevant solution to be analysed. The test based on inhibition of bacterial bioluminescence is used to determine its potential toxicity. It can be diluted to assess its dose-response relationship.

123  
 124  
 125  
 126  
 127  
 128  
 129  
 130

131 **2. Software description**

132

133 This R-Shiny package is written in the R [19] and Shiny framework [20] and hosted on a github repository  
134 (available at [https://github.com/bbellier/biolutoxR\\_package](https://github.com/bbellier/biolutoxR_package)). This Shiny App used only a few  
135 dependencies (Table 2). To use this Shiny application package, the following commands can be used:

136

137 (1) To install the *biolutoxR* package:

138

139 *R> install.packages("devtools")*

140 *R> devtools::install\_github("bbellier/biolutoxR\_package")*

141 *R> library(biolutoxR)*

142

143 (2) To run the *biolutoxR* Shiny App:

144

145 *R> run.biolutoxR()*

146

147 (3) To run the *biolutoxR* Shiny App example:

148

149 *R> example.biolutoxR()*

150

151 (4) To cite the *biolutoxR* package:

152

153 *R> citation.biolutoxR()*

154

155

156 **Table 2.** *biolutoxR* Shiny App dependencies.

157

Dependencies	Functions	CRAN Reference
bslib	Use to customize Shiny CSS styling	<a href="https://cran.r-project.org/web/packages/bslib/index.html">https://cran.r-project.org/web/packages/bslib/index.html</a>
DT	Use to table management in Shiny	<a href="https://cran.r-project.org/web/packages/DT/index.html">https://cran.r-project.org/web/packages/DT/index.html</a>
openxlsx	Use to generate and read .xlsx file	<a href="https://cran.r-project.org/web/packages/openxlsx/index.html">https://cran.r-project.org/web/packages/openxlsx/index.html</a>
shiny	Use to generate an R-Shiny application	<a href="https://cran.r-project.org/web/packages/shiny.exe/index.html">https://cran.r-project.org/web/packages/shiny.exe/index.html</a>
shinythemes	Use to customize Shiny aspect	<a href="https://cran.r-project.org/web/packages/shinythemes/index.html">https://cran.r-project.org/web/packages/shinythemes/index.html</a>
shinyMatrix	Use to generate matrix in Shiny	<a href="https://cran.r-project.org/web/packages/shinyMatrix/index.html">https://cran.r-project.org/web/packages/shinyMatrix/index.html</a>
tidyverse	Useful set of packages (e.g. ggplot package for graphical data visualization)	<a href="https://cran.r-project.org/web/packages/tidyverse/index.html">https://cran.r-project.org/web/packages/tidyverse/index.html</a>

158

159 The operating principle of the *biolutoxR* package application is summarized in Figure 2. Therefore, this  
160 Shiny App is divided into different tabs that can be navigated using a navigation bar (Fig. 3). The application  
161 is therefore composed of 8 tabs: “*Presentation*”, “*Scheme*”, “*Time 5*”, “*Time 15*”, “*Time 30*”, “*Table*”,  
162 “*Plot*”, “*Exit*” and 4 main uses: R-Shiny package presentation, bacterial bioluminescence toxicity test data  
163 input, final data table construction and final data visualization.

164  
165 The “*Presentation*” tab is necessary to help users understand (1) how the bacterial bioluminescence toxicity  
166 test works and (2) the relevance and operation of the Shiny application. Indeed, the first paragraph  
167 succinctly describes the main stages of a bacterial bioluminescence toxicity test, its usefulness, and the  
168 analysis that can be made of it. The second paragraph guides the user through the mechanisms and steps in  
169 inputting data and understanding the final constructed variables.

170  
171 The ‘*Scheme*’ tab contains a matrix used to define the identifiers of the manipulation (Fig. 4). In this tab,  
172 the user must define the name of the manipulation (called ‘*manip\_name*’ in the final data table), the number  
173 of matrices required (i.e. one matrix represents one set, ‘*Number of matrices*’), the size of the matrices  
174 (‘*Number of rows*’ and ‘*Number of columns*’, which allow to define the position of a sample in the matrix,  
175 which is given in the ‘*id*’ variable in the final data table) and the name of the negative control (‘*Name of*  
176 *negative control*’), which is then used to correct the bioluminescence values (called *biolu* in the ‘*Time*’ tabs  
177 and in the final data table) in the final data table. Each cell is filled to position and characterise the samples  
178 in the matrix. Each cell of the matrix must contain values in this format: ‘*sol,dil,rep\_bio,sol\_sh*’; where  
179 ‘*sol*’ corresponds to the full name of the solution (e.g. tire and road wear particles of France), ‘*dil*’  
180 corresponds to the dilution level (from 0% to 100%), ‘*rep\_bio*’ corresponds to the biological replicate and  
181 ‘*sol\_sh*’ corresponds to the short name of the solution (e.g. TRWP-FR).


182  
183 The ‘*Time 1*’, ‘*Time 2*’ and ‘*Time 3*’ tabs each contain a matrix which works in the same way (Fig. 5). For  
184 each *Time* tab, the exposure time (‘*Time (in min)*’) is entered by the user to be considered in the final data  
185 table. Each cell in these matrices is filled with bioluminescence data from the bioluminescence bacterial  
186 toxicity test.

187  
188 The ‘*Table*’ tab displays the final data table obtained, which can be exported in ‘.csv’ or ‘.xlsx’ format. All  
189 the data generated by completing the ‘*Scheme*’ and ‘*Time*’ matrices is included in this final data table:  
190 ‘*manip\_name*’, ‘*set*’, ‘*time*’, ‘*id*’, ‘*sol*’, ‘*sol\_sh*’, ‘*dil*’, ‘*rep\_bio*’, ‘*biolu*’, ‘*biolu\_mean\_neg\_control*’,  
191 ‘*biolu\_corr*’, ‘*perc\_inhib*’, ‘*perc\_inhib\_corr*’). The final value used to display the results is  
192 ‘*perc\_inhib\_corr*’ (i.e. corrected inhibition percentage, i.e. bioluminescence inhibition corrected values),  
193 which requires the calculation of ‘*biolu*’ (i.e. bioluminescence measured values),  
194 ‘*biolu\_mean\_neg\_control*’ (i.e. mean of bioluminescence measured values for the negative control sample),  
195 ‘*biolu\_corr*’ (i.e. bioluminescence corrected values) and ‘*perc\_inhib*’ (i.e. bioluminescence inhibition  
196 calculated values).


197  
198 The ‘*Plot*’ tab provides a simple and dynamic way of viewing and modifying plots and generating  
199 standardised toxicity plots/data (i.e. the dose-response curve and EC<sub>x</sub>). Firstly, in ‘*Final data plot*’ part, the  
200 user can modify the visualisation of the results as desired (choice of x variable, y variable, time, solutions,  
201 dilutions and colour variables). Secondly, in ‘*Dose-response curve*’ part, the dose-response curve is plotted  
202 (by setting the the *c* parameter to 0 and *d* parameter to 100), using the appropriate dilution(s) and substance

203 of interest. Furthermore, the equation of the dose-response curve, the significance of the parameters (b, i.e.  
 204 the curvature parameter and e, i.e. the  $EC_{50}$ ) and the Pearson's linear coefficient of determination ( $R^2$ , i.e.  
 205 measure to assess the efficiency of a linear regression model; an  $R^2$  of 1 indicates a robust model) are  
 206 calculated. Also, any desired  $EC_X$  (with its standard deviation and 95% confidence interval) could be  
 207 calculated and printed.

208  
 209 The “Exit” tab serves to quit the Shiny App.  
 210  
 211



212  
 213 **Figure 2.** biolutoxR package structure with its main functionalities.  
 214  
 215



216  
 217 **Figure 3.** biolutoxR Shiny app navigation bar.  
 218  
 219

	1	2	3	4	5
A	sol,dil,rep_bio,sol_sh	sol,dil,rep_bio,sol_sh	sol,dil,rep_bio,sol_sh	sol,dil,rep_bio,sol_sh	sol,dil,rep_bio,sol_sh
B	sol,dil,rep_bio,sol_sh	sol,dil,rep_bio,sol_sh	sol,dil,rep_bio,sol_sh	sol,dil,rep_bio,sol_sh	sol,dil,rep_bio,sol_sh
C	sol,dil,rep_bio,sol_sh	sol,dil,rep_bio,sol_sh	sol,dil,rep_bio,sol_sh	sol,dil,rep_bio,sol_sh	sol,dil,rep_bio,sol_sh
D	sol,dil,rep_bio,sol_sh	sol,dil,rep_bio,sol_sh	sol,dil,rep_bio,sol_sh	sol,dil,rep_bio,sol_sh	sol,dil,rep_bio,sol_sh
E	sol,dil,rep_bio,sol_sh	sol,dil,rep_bio,sol_sh	sol,dil,rep_bio,sol_sh	sol,dil,rep_bio,sol_sh	sol,dil,rep_bio,sol_sh
F	sol,dil,rep_bio,sol_sh	sol,dil,rep_bio,sol_sh	sol,dil,rep_bio,sol_sh	sol,dil,rep_bio,sol_sh	sol,dil,rep_bio,sol_sh

220  
 221 **Figure 4.** biolutoxR app Shiny “Scheme” tab for an 5x6 matrix.  
 222

	1	2	3	4	5
A	biolu	biolu	biolu	biolu	biolu
B	biolu	biolu	biolu	biolu	biolu
C	biolu	biolu	biolu	biolu	biolu
D	biolu	biolu	biolu	biolu	biolu
E	biolu	biolu	biolu	biolu	biolu
F	biolu	biolu	biolu	biolu	biolu

223  
 224 **Figure 5.** biolutoxR app Shiny “Time” tab for an 5x6 matrix.

225  
226  
227  
228  
229  
230  
231  
232  
233  
234  
235  
236  
237  
238  
239  
240  
241  
242  
243  
244  
245  
246  
247  
248  
249  
250  
251

### 3. Illustrative example

The example used is not based on actual data. This example can be found in the function: *example.biolutoxR()*.

In a geographical area, new agricultural and industrial practices were introduced close to a river (Fig. 6). Ecological disasters followed rapidly afterwards, with fish species seeing their populations collapse.

Thus, local authorities have called on the scientific community to study and understand this phenomenon. To understand the toxic effect of new molecules used in agriculture and industry close to the river, scientists proposed the use of a bioassay: the toxicity test based on bacterial bioluminescence inhibition; using samples of the new agrochemicals and industrial effluent at various relevant concentrations (Fig. 6).

Following the test, the bacterial bioluminescence values were recovered and entered into the Shiny application offered by the *biolutoxR* package. A cleaned final table was exported with the test results (Fig. 7) and a graph was obtained to analyse the global situation of this locality (Fig. 8). The graph showed that even at low dilution, the percentage inhibition of bioluminescence indicates high toxicity for both effluents. The scientists then focused on substance B (i.e. industrial effluent), which appears to be more toxic, to obtain its  $EC_{50}$  value from its dose-response curve. An  $EC_{50}$  of 12.33% (with a standard deviation of 0.86 and a 95% confidence interval of 9.58-15.08) was obtained for the industrial effluent. In other words, the concentration required to obtain a 50% response in exposed organisms compared with the negative control was 12.33% (Fig. 9). Furthermore, the value of the  $R^2$  (0.99) and the significance of the  $EC_{50}$  parameter ( $< 0.05$ ) can be verified to ensure the robustness of the  $EC_x$  value obtained.




Figure 6. Context of the example study.

252  
253  
254  
255  
256



**Final data table:**

Download CSV Download XLSX

Show 30 entries Search:

	manip_name	set	time	id	sol	sol_sh	dil	rep_bio	biolu	biolu_mean_neg_control	biolu_corr	perc_inhib	perc_inhib_corr
1	Example	Set 1	5	A1	NegControl	Neg	100	1	100	96.2	103.95	-3.95	0
2	Example	Set 2	5	A1	NegControl	Neg	100	2	100	96.2	103.95	-3.95	0
3	Example	Set 1	5	A2	NegControl	Neg	100	1	96	96.2	99.79	0.21	0.21
4	Example	Set 2	5	A2	NegControl	Neg	100	2	96	96.2	99.79	0.21	0.21
5	Example	Set 1	5	A3	NegControl	Neg	100	1	100	96.2	103.95	-3.95	0
6	Example	Set 2	5	A3	NegControl	Neg	100	2	100	96.2	103.95	-3.95	0
7	Example	Set 1	5	A4	NegControl	Neg	100	1	95	96.2	98.75	1.25	1.25
8	Example	Set 2	5	A4	NegControl	Neg	100	2	95	96.2	98.75	1.25	1.25
9	Example	Set 1	5	A5	NegControl	Neg	100	1	90	96.2	93.56	6.44	6.44

Showing 1 to 30 of 90 entries Previous 1 2 3 Next

Figure 7. Screenshot of the final data table.

**Final data plot:**




Figure 8. Screenshot of the simple and dynamic plot proposed by the application.

**Dose-Response curve:**




Figure 9. Screenshot of the dose-response curve and EC<sub>x</sub> results proposed by the application.

257  
258  
259

260  
261  
262

263  
264  
265  
266  
267

Further information is available on the website dedicated to the *biolutoxR* package (available at [https://bbellier.github.io/biolutoxR\\_website/](https://bbellier.github.io/biolutoxR_website/)).

268  
269  
270  
271  
272  
273  
274  
275  
276  
277  
278  
279  
280  
281  
282  
283  
284  
285  
286  
287  
288  
289  
290  
291  
292  
293  
294  
295

#### 4. Impact

The *biolutoxR* package simplify data analysis of toxicity tests based on bacterial bioluminescence inhibition for beginners and experimented users. Indeed, an introduction to this type of test is included. Furthermore, this package providing an automated, simple and dynamic tool to enter, recover and visualize results. In addition to providing an example, an introduction to the application enables users to quickly understand its functioning. Widespread use of this tool, particularly through the construction of a dose-response curve and  $EC_x$  calculations, will reduce the time allocated to data cleaning and calculations, increase the reproducibility and comparability of data and facilitate access to ecotoxicological data.

The *biolutoxR* package can be used for all environmental compartments (air, sediment, freshwater and marine water) and applications can vary from the assessment of the toxicity of a single molecule (e.g. N-(1,3-dimethylbutyl)-N'-phenyl-p-phenylenediamine), of a mixture of molecules (e.g. tire and road wear particles) or of a matrix (e.g. stormwater runoff). In this way, this tool can contribute to improving our knowledge of toxicity mechanisms, particularly by supplementing the ecotoxicity databases available with toxicity values ( $EC_x$ ), which are currently still very limited in terms of quality and quantity [21].

#### 5. Conclusion

Overall, the *biolutoxR* package facilitates data analysis (data cleaning, results plotting, toxicity data accessing) while retaining the use of the classic matrix format for a toxicity test based on bacterial bioluminescence inhibition. More concretely, the advantages of this toxicity test (speed, reproducibility and low cost) ensure the relevance of this tool for improving ecotoxicological knowledge for a wide range of applications. To our knowledge, this R-Shiny application represents the first initiative to generalize flexible data analysis for these bioassays, improving the reproducibility and comparability of the data recovered following these tests.

296 **Bibliography.**  
297

- 298 [1] Richardson K, Steffen W, Lucht W, Bendtsen J, Cornell SE, Donges JF, et al. Earth beyond six of  
299 nine planetary boundaries. *Sci Adv* 2023;9:eadh2458. <https://doi.org/10.1126/sciadv.adh2458>.
- 300 [2] Ali H, Khan E. Environmental chemistry in the twenty-first century. *Environ Chem Lett*  
301 2017;15:329–46. <https://doi.org/10.1007/s10311-016-0601-3>.
- 302 [3] Borgwardt F, Robinson L, Trauner D, Teixeira H, Nogueira AJA, Lillebø AI, et al. Exploring  
303 variability in environmental impact risk from human activities across aquatic ecosystems. *Science of*  
304 *The Total Environment* 2019;652:1396–408. <https://doi.org/10.1016/j.scitotenv.2018.10.339>.
- 305 [4] Schuijt LM, Peng F-J, van den Berg SJP, Dingemans MML, Van den Brink PJ. (Eco)toxicological  
306 tests for assessing impacts of chemical stress to aquatic ecosystems: Facts, challenges, and future.  
307 *Science of The Total Environment* 2021;795:148776.  
308 <https://doi.org/10.1016/j.scitotenv.2021.148776>.
- 309 [5] Azizullah A, Häder D-P. A comparison of commonly used and commercially available bioassays for  
310 aquatic ecosystems. *Bioassays*, Elsevier; 2018, p. 347–68. [https://doi.org/10.1016/B978-0-12-](https://doi.org/10.1016/B978-0-12-811861-0.00017-6)  
311 [811861-0.00017-6](https://doi.org/10.1016/B978-0-12-811861-0.00017-6).
- 312 [6] Ma XY, Wang XC, Ngo HH, Guo W, Wu MN, Wang N. Bioassay based luminescent bacteria:  
313 Interferences, improvements, and applications. *Science of The Total Environment* 2014;468–469:1–  
314 11. <https://doi.org/10.1016/j.scitotenv.2013.08.028>.
- 315 [7] Neale PA, O'Brien JW, Glauch L, König M, Krauss M, Mueller JF, et al. Wastewater treatment  
316 efficacy evaluated with in vitro bioassays. *Water Research X* 2020;9:100072.  
317 <https://doi.org/10.1016/j.wroa.2020.100072>.
- 318 [8] Masner P, Javůrková B, Bláha L. Rapid in situ toxicity testing with luminescent bacteria  
319 *Photorhabdus luminescens* and *Vibrio fischeri* adapted to a small portable luminometer. *Environ Sci*  
320 *Pollut Res* 2017;24:3748–58. <https://doi.org/10.1007/s11356-016-8096-9>.
- 321 [9] Parvez S, Venkataraman C, Mukherji S. A review on advantages of implementing luminescence  
322 inhibition test (*Vibrio fischeri*) for acute toxicity prediction of chemicals. *Environment International*  
323 2006;32:265–8. <https://doi.org/10.1016/j.envint.2005.08.022>.
- 324 [10] Abbas M, Adil M, Ehtisham-ul-Haque S, Munir B, Yameen M, Ghaffar A, et al. *Vibrio fischeri*  
325 bioluminescence inhibition assay for ecotoxicity assessment: A review. *Science of The Total*  
326 *Environment* 2018;626:1295–309. <https://doi.org/10.1016/j.scitotenv.2018.01.066>.
- 327 [11] Backhaus T, Grimme LH. The toxicity of antibiotic agents to the luminescent bacterium *Vibrio*  
328 *fischeri*. *Chemosphere* 1999;38:3291–301. [https://doi.org/10.1016/S0045-6535\(98\)00560-8](https://doi.org/10.1016/S0045-6535(98)00560-8).
- 329 [12] Muneeswaran T, Kalyanaraman N, Vennila T, Rajesh Kannan M, Ramakritinan CM. Rapid  
330 assessment of heavy metal toxicity using bioluminescent bacteria *Photobacterium leiognathi* strain  
331 GoMGm1. *Environ Monit Assess* 2021;193:109. <https://doi.org/10.1007/s10661-021-08860-2>.
- 332 [13] Thomulka KW, Schroeder JA, Lange JH. Use of *Vibrio harveyi* in an aquatic bioluminescent toxicity  
333 test to assess the effects of metal toxicity: Treatment of sand and water-buffer, with and without  
334 EDTA. *Environ Toxicol Water Qual* 1997;12:343–8. [https://doi.org/10.1002/\(SICI\)1098-](https://doi.org/10.1002/(SICI)1098-2256(1997)12:4<343::AID-TOX8>3.0.CO;2-8)  
335 [2256\(1997\)12:4<343::AID-TOX8>3.0.CO;2-8](https://doi.org/10.1002/(SICI)1098-2256(1997)12:4<343::AID-TOX8>3.0.CO;2-8).
- 336 [14] Girotti S, Ferri EN, Fumo MG, Maiolini E. Monitoring of environmental pollutants by bioluminescent  
337 bacteria. *Analytica Chimica Acta* 2008;608:2–29. <https://doi.org/10.1016/j.aca.2007.12.008>.
- 338 [15] Si H, Zhou G, Luo Y, Wang Z, Pan X, Dao G. Hormesis in the Assessment of Toxicity Assessment  
339 by Luminescent Bacterial Methods. *Toxics* 2024;12:596. <https://doi.org/10.3390/toxics12080596>.

- 340 [16] Westlund P, Nasuhoglu D, Isazadeh S, Yargeau V. Investigation of Acute and Chronic Toxicity  
341 Trends of Pesticides Using High-Throughput Bioluminescence Assay Based on the Test Organism  
342 *Vibrio fischeri*. Arch Environ Contam Toxicol 2018;74:557–67. [https://doi.org/10.1007/s00244-017-](https://doi.org/10.1007/s00244-017-0483-9)  
343 0483-9.
- 344 [17] Tothill IE, Turner APF. Developments in bioassay methods for toxicity testing in water treatment.  
345 TrAC Trends in Analytical Chemistry 1996;15:178–88. [https://doi.org/10.1016/0165-](https://doi.org/10.1016/0165-9936(96)80640-6)  
346 9936(96)80640-6.
- 347 [18] Kokkali V, Van Delft W. Overview of commercially available bioassays for assessing chemical  
348 toxicity in aqueous samples. TrAC Trends in Analytical Chemistry 2014;61:133–55.  
349 <https://doi.org/10.1016/j.trac.2014.08.001>.
- 350 [19] R Core Team. R: A language and environment for statistical computing. Vienne, Austria: R  
351 foundation for statistical computing; 2021.
- 352 [20] Chang W, Cheng J, Allaire JJ, Sievert C, Schloerke B, Xie Y, et al. Shiny: Web Application  
353 Framework for R. 2024.
- 354 [21] Bellier B, Bancel S, Rochard É, Cachot J, Geffard O, Villeneuve B. Assessment of the impact of  
355 chemical pollution on endangered migratory fish in two major rivers of France, including spawning  
356 grounds. Science of The Total Environment 2024;931:172748.  
357 <https://doi.org/10.1016/j.scitotenv.2024.172748>.  
358

359 **Author contact.** Coralie Le Picard; E-mail address: coralie.lepicard@yahoo.com; ORCID: 0009-0000-  
360 5928-7681.

361  
362 **Credit authorship contribution statement:** *Coralie Le Picard*: Conceptualization, Writing & Editing.  
363 *Benjamin Bellier*: Conceptualization, Writing & Editing.

364  
365 **Declaration of competing interest.** The authors declare that they have no known competing financial  
366 interests or personal relationships that could have appeared to influence the work reported in this paper.

367  
368 **Data availability.** The package is accessible on the github of Bellier Benjamin  
369 ([https://github.com/bbellier/biolutoxR\\_package](https://github.com/bbellier/biolutoxR_package)) and on the website dedicated to the *biolutoxR* package  
370 ([https://bbellier.github.io/biolutoxR\\_website/](https://bbellier.github.io/biolutoxR_website/)).

371  
372 **Acknowledgments.** This work was funded by the PLASTyre project funded by a partnership between:  
373 “Région Nouvelle Aquitaine, Département de la Charente Maritime”, “Agence de l’Eau Adour-Garonne”  
374 and “Qualyse”. The funding was also supported in part by: “La Rochelle University”, “University of  
375 Bordeaux”, “Centre National de la Recherche Scientifique”. We thank anonymous reviewers for their kind  
376 and insightful comments on a previous version of the manuscript and package, Quentin Petitjean for helpful  
377 discussions at various stages of this package development.

378  
379

380 **Supplementary materials.**


381

382 *Supplementary materials n°1.*

383

384 General reaction equation of bacterial bioluminescence:

385



387

388 In other terms:

389

- 390 • Long-chain aldehydes + reduced flavin mononucleotide + molecular oxygen — luciferase → long-
- 391 chain acids + oxidized flavin mononucleotide + water + light emission

392

393